

VELA[®] ventilator

Mode

Visual display of current mode with easy and flexible access to all modes of operation

User selectable patient monitor values

Five monitored parameters selectable from 19 options are displayed continuously on main screen

Ventilation settings


Ventilator settings are displayed on all screens

Main screen

An active matrix, full color LCD touch screen features scaleable graphic displays

Nebulizer

Synchronized on inspiration, nebulization time can be set in increments of one minute for up to a maximum of 60 minutes



Main
Access additional menus easily from the main screen

Alarm silence
1 minute

Alarm message window
Alarm messages are displayed visually with audible sound when active according to priority

Manual breath
100% O₂
3 minutes

Intuitive operation
The Simple Touch™ user interface in the VELA[®] ventilator reduces training time, produces consistency in method and enhances clinician confidence.

WARNING—U.S. Federal Law restricts this device to sale by or on the order of a physician.



CareFusion
Yorba Linda, CA

carefusion.com

CareFusion
22745 Savi Ranch Parkway
Yorba Linda, CA 92887
800.231.2466 toll-free
714.283.2228 tel
714.283.8493 fax

CareFusion Germany GmbH 234
Leibnizstrasse 7
D-97204 Hoechberg
Germany
+49.0.931.4972.0 tel
+49.0.931.4972.423 fax

